

Anti-NSL1 antibody (1-281) (STJ29169)

STJ29169

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/ICC:1:50-1:100
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25936
Gene Symbol	NSL1
Uniprot ID	NSL1_HUMAN
Immunogen	
Immunogen Region	1-281
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-281 of human NSL1 (NP_056286.3).
Immunogen Sequence	MAGSPELVLDPPWDKELAA GTESQALVSATPREDFRVRC TSKRAVTEMLQLCGRFVQKL GDALPEEIREPALRDAQWTF ESAVQENISINGQAWQEASD NCFMDSDIKVLEDQFDEIIV DIATKRKQYPRKILECVIKT IKAKQEILKQYHPVVHPLDL KYDPDPAPHMENLKCRGETV AKEISEAMKSLPALIEQGEG FSQVLRMQPVIHLQRIHQEV FSSCHRKPDAPENFITQI



Immunofluorescence analysis of A549 cells using NSL1 antibody (STJ29169).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081