

Anti-SH3GLB1 antibody (1-365) (STJ29184)

STJ29184

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51100
Gene Symbol	SH3GLB1
Uniprot ID	SHLB1_HUMAN
Immunogen	
Immunogen Region	1-365
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-365 of human Bif-1 (NP_057093.1).
Immunogen Sequence	MNIMDFNVKKLAADAGTFLS RAVQFTEKLGQAETELDA HLENLLSKAECTKIWTEKIM KQTEVLLQPNPNARIEEFVY EKLDKAPSRINPELLGQY MIDAGTEFGPGTAYGNALIK CGETQKRIGTADRELIQTS A LNFLTPLRNFIEGDYKTIK ERKLLQNKRLDLDAKTRLK KAKAAETRNSSEQLRITQS EFDRQAEITRLLLEGISSTH AHHLRCLNDFVEAQMTYYA



Immunofluorescence analysis of A549 cells using Bif-1 Rabbit polyclonal antibody (STJ29184). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.