

Anti-ANAPC5 antibody (1-250) (STJ29189)

STJ29189

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51433
Gene Symbol	ANAPC5
Uniprot ID	APC5_HUMAN
Immunogen	
Region	1-250
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ANAPC5 (NP_057321.2).
Immunogen Sequence	MASVHESLYFNPMMTNGVVH ANVFGIKDWTPYKIAVLVL LNEMSRGTGEGAVSLMERRRL NQLLLPLLQGPDTLSKLYK LIEESCPLANSVQIRIKLM AEGELKDMEQFFDDLSDSFS GTEPEVHKTSVVGFLRLRHM LAYSKLFSQVFKLYTALQQ YFQNGEKKTVEDADMELTSR DEGERKMEKEELDVSVREEE VSCSGPLSQKQAEFFLSQQA SLLKNDETALTPASLQKE



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081