

Anti-PINK1 antibody (282-581) (STJ29211)

STJ29211

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65018
Gene Symbol	PINK1
Uniprot ID	PINK1_HUMAN
Immunogen	
Immunogen Region	282-581
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 282-581 of human PINK1 (NP_115785.1).
Immunogen Sequence	TSSVPLLPGALVDYDPDLPS RLHPEGLGHGRTLFLVMKNY PCTLRQYLCVNTSPRLAAM MLLQLLEGVDHLVQQGIAHR DLKSDNILVELDPDGCPWLVIADFGCCLADESIGLQLPFS SWYVDRGGNGCLMAPEVSTA RPGPRAVIDYSKADAWAVGA IAYEIFGLVNPFGQGKAHL ESRSYQEAQLPALPESVPPD VRQLVRALLQREASKRPSAR VAANVLHLSLWGEHILALK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081