

Anti-LELP1 antibody (1-98) (STJ29228)

STJ29228

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	149018
Gene Symbol	LELP1
Uniprot ID	LELP1_HUMAN
Immunogen	
Immunogen Region	1-98
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-98 of human LELP1 (NP_001010857.1).
Immunogen Sequence	MSSDDKSKSNDPKTEPKNCD PKCEQKCESKCQPSCLKLL QRCFEKCPWEKCPAPPKCLP CPSQSPSSCPPQPCTKPCPP KCPSSCPHACPPPCPPPE