

Anti-COMMD1 antibody (1-190) (STJ29229)

STJ29229

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	150684
Gene Symbol	COMMD1
Uniprot ID	COMD1_HUMAN
Immunogen	
Immunogen Region	1-190
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human COMMD1 (NP_689729.1).
Immunogen Sequence	MAAGELEGGKPLSGLLNALA QDTFHGYPGITEELLRSQLY PEVPPPEEFRPFLAKMRGILK SIASADMDFNQLAFLTAQT KKQGGITSDQAAVISKFWKS HKTKIRESLMNQSRWNSGLR GLSWRVGKSGSRHSAQIHT PVAIIELGKYGGQESEFLC LEFDEVKNQILKTLSEVEE SISTLISQPN



Immunofluorescence analysis of U2OS cells using COMMD1 antibody (STJ29229) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of L929 cells using COMMD1 antibody (STJ29229) at dilution of 1:100. Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081