

# Anti-PHIP antibody (1542-1821) (STJ29287)

STJ29287

## GENERAL INFORMATION


|                          |                       |
|--------------------------|-----------------------|
| <b>Product Type</b>      | Primary antibodies    |
| <b>Short Description</b> |                       |
| <b>Applications</b>      | WB/IHC-P/IF/ICC/ELISA |
| <b>Host/Source</b>       | Rabbit                |
| <b>Reactivity</b>        | Human/Mouse           |

## PRODUCT PROPERTIES


|                      |                                                                                                                                                                                 |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>     | Polyclonal                                                                                                                                                                      |
| <b>Clone ID</b>      |                                                                                                                                                                                 |
| <b>Concentration</b> | Lot specific                                                                                                                                                                    |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                    |
| <b>Purification</b>  | Affinity purification                                                                                                                                                           |
| <b>Dilution</b>      | WB:1:500-1:2000                                                                                                                                                                 |
| <b>Range</b>         | IHC-P:1:50-1:200<br>IF/ICC:1:50-1:100<br>ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>   | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.                                                                                                                              |
| <b>Isotype</b>       | IgG                                                                                                                                                                             |
| <b>Storage</b>       | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                                                  |
| <b>Instruction</b>   |                                                                                                                                                                                 |

## TARGET INFORMATION


|                           |                                                                                                                                                                                                                                                               |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 55023                                                                                                                                                                                                                                                         |
| <b>Gene Symbol</b>        | PHIP                                                                                                                                                                                                                                                          |
| <b>Uniprot ID</b>         | PHIP_HUMAN                                                                                                                                                                                                                                                    |
| <b>Immunogen</b>          |                                                                                                                                                                                                                                                               |
| <b>Region</b>             | 1542-1821                                                                                                                                                                                                                                                     |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 1542-1821 of human PHIP (NP_060404.3).                                                                                                                                          |
| <b>Immunogen Sequence</b> | KTILENSVKHISKALNTLSPP GQSSFSHGTRNNSAKENMEK EKPVKRKMKSSVLPKASTLS KSSAVIEQGDCCKNNALVPGT<br>IQVNGHGGQPSKLVKRGPGR KPKVEVTNSGEIHKRGR KPKKLQYAKPEDLEQNNVHP IRDEVLPSTCNFLSETNNV<br>KEDLLQKKNRGGKPKRKMK TQKLDADLLVPASVKVLRSS NRKIDDDPDEEEEFEEELKG SEPHMRTRNQGRRTAFYNE |




Immunofluorescence analysis of MCF-7 cells using PHIP antibody (STJ29287).



Immunohistochemistry analysis of paraffin-embedded mouse tumor using PHIP antibody (STJ29287) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human kidney cancer using PHIP antibody (STJ29287) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human well-differentiated squamous skin carcinoma using PHIP antibody (STJ29287) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081