

Anti-TMLHE antibody (167-376) (STJ29451)

STJ29451

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55217
Gene Symbol	TMLHE
Uniprot ID	TMLH_HUMAN
Immunogen	
Immunogen Region	167-376
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 167-376 of human TMLHE (NP_001171726.1).
Immunogen Sequence	CQSFLLETNEGLKKFLQNFL YGIAFVENVPPTQEHTKLA ERISLI RETIYGRMWYFTSD FSRGDTAYTKLALDRHTD TT YFQEPGCIQVFHCLKHEGTG GRTLLVDGFYAAEQVLQKAP EEFELLSKVPLKHEYIEDVG ECHNHMIGIGPVLNIYPWNK ELYLIRLFKEKQNTVNRQWN SSLQCDIPERILTYRHFVSG TSIEHRGSLI



Immunofluorescence analysis of MCF-7 cells using TMLHE antibody (STJ29451). Blue: DAPI for nuclear staining.