

Anti-BMP15 antibody (143-392) (STJ29460)
STJ29460

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:100 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9210
Gene Symbol	BMP15
Uniprot ID	BMP15_HUMAN
Immunogen	
Immunogen Region	143-392
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 143-392 of human BMP15 (NP_005439.2).
Immunogen Sequence	LTRFNLSCHVWVQKNPTN HFPSSGDSKPSLMSNAWK EMDITQLVQQRFWNNKGHRI LRLRFMCQQKQDSGGLELWH GTSSLDIAFLLLYFNDTHKS IRKAKFLPRGMEEFMERESL LRRTRQADGISA EVTASSSK HSGPENNQCSLHPFQISFRQ LGWDHWIAPPFYTPNYCKG TCLRVLRDLNLSNPHAIQN LINQLVDQSVPRPSCVPYKY VPISVLMIEANGSILYKEY



Immunofluorescence analysis of MCF7 cells using BMP15 Rabbit polyclonal antibody (STJ29460). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using BMP15 Rabbit polyclonal antibody (STJ29460) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081