

Anti-ATG9B antibody (715-924) (STJ29543)

STJ29543

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	285973
Gene Symbol	ATG9B
Uniprot ID	ATG9B_HUMAN
Immunogen	
Immunogen Region	715-924
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 715-924 of human ATG9B (NP_001303985.1).
Immunogen Sequence	HPLWRPPGHSSKFLGHLWGR VQQDAAAWGATSARGPSTPG VLSNCTSPLEAFLANLFVH PLLPPRDLSP TAPCPAAATA SLLASISRIAQDPSSVSPGG TGGQKLAQLPELASAEMSLH VIYLHQLHQQQQQEPWGEA AASILSRPCSSSPQPPSPDE EKPSWSSDGGSPASSPRQWV GTQKARNLFPGGFQVTTDTQ KEPDRASCTD



Immunofluorescence analysis of human placenta cells using ATG9B rabbit polyclonal antibody (STJ29543) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081