

Anti-TEAD3 antibody (130-330) (STJ29590)

STJ29590

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	7005
Gene Symbol	TEAD3
Uniprot ID	TEAD3_HUMAN
Immunogen	
Region	130-330
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 130-330 of human TEAD3 (NP_003205.2).
Immunogen Sequence	SSAQIVSASVLQNKFSPPSP LPQAVFSTSSRFWSSPPLLG QQPGPSQDIKPFAPAYPIQ PPLPPTLSSYEPLAPLPSAA ASVPVWQDRTIASSRLRLLE YSAFMEVQRDPDTYSKHLFV HIGQTNPAFSDPPEAVDVR QIYDKFPEKKGGKLKELYEK PPNAFLVKFWADLNSTIQE GPQAFYGVSSQYSSADSM TI S




Confocal immunofluorescence analysis of U-2 OS cells using TEAD3 Rabbit polyclonal antibody (STJ29590) at dilution of 1:400. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of A549 cells using TEAD3 Rabbit polyclonal antibody (STJ29590). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution.



Immunohistochemistry analysis of TEAD3 in paraffin-embedded human placenta using TEAD3 Rabbit polyclonal antibody (STJ29590) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using TEAD3 Rabbit polyclonal antibody (STJ29590) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ29590) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081