

Anti-TNFRSF4 antibody (29-214) (STJ29635)
STJ29635

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7293
Gene Symbol	TNFRSF4
Uniprot ID	TNR4_HUMAN
Immunogen	
Immunogen Region	29-214
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 29-214 of human CD134/OX40 (NP_003318.1).
Immunogen Sequence	LHCVGDTPSNDRCHECRP GNGMVSRCRSQNTVCRPCG PGFYNDVSSKPCKPCTWCN LRSGSERKQLCTATQDTVCR CRAGTQPLDSYKPGVDCAPC PPGHFSPGDNQACKPWTNCT LAGKHTLQPASNSSDAICED RDPPATQPQETQGPPARPIT VQPTEAWPRTSQGPSTRPVE VPGGRA

