

Anti-FANCM antibody (390-660) (STJ29739)

STJ29739

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:100
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57697
Gene Symbol	FANCM
Uniprot ID	FANCM_HUMAN
Immunogen	
Immunogen Region	390-660
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 390-660 of human FANCM (NP_065988.1).
Immunogen Sequence	MDGKGMTRSKNELGRNEDF MKLYNHLECMFARTRTSAN GISAIQQGDKNKKFVYSHPK LKKEEVIEHFKSWNAENT TEKKRDETRVMIFSSFRDSV QEIAEMLSQHQPIIRVMTFV GHASGKSTKGFTQKEQLEV KQFRDGGYNTLVSTCVGEEG LDIGEVDLIICFDSQKSPIR LVQRMGRTGRKRQGRIVIL SEGREERIYNQSQSNKRSIY KAISNRQVLHIFYQRSPRM



Immunofluorescence analysis of A549 cells using FANCM antibody (STJ29739).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081