

## Anti-Albumin antibody [6502] (STJ400008)

STJ400008

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Albumin is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	6502
<b>Concentration</b>	2.5 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purity accurate to greater than or equal to 95%
<b>Dilution Range</b>	
<b>Formulation</b>	0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>Isotype</b>	IgG2a
<b>Storage</b>	36 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21
<b>Instruction</b>	days OK, but some charge alterations +35 °C, 7 days OK, but some charge alterations +35 °C, 21 days OK, but some charge alter

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">213</a>
<b>Gene Symbol</b>	<a href="#">ALB</a>
<b>Uniprot ID</b>	<a href="#">ALBU_HUMAN</a>
<b>Immunogen</b>	human albumin
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	