

**Anti-C-Peptide antibody [9101] (STJ400017)**

STJ400017

**GENERAL INFORMATION**

|                          |                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                           |
| <b>Short Description</b> | Mouse monoclonal antibody anti-C-Peptide is suitable for use in ELISA research applications. |
| <b>Applications</b>      | ELISA                                                                                        |
| <b>Host/Source</b>       | Mouse                                                                                        |
| <b>Reactivity</b>        | Human                                                                                        |

**PRODUCT PROPERTIES**

|                       |                                                                                                                              |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Monoclonal                                                                                                                   |
| <b>Clone ID</b>       | 9101                                                                                                                         |
| <b>Concentration</b>  | 5.0 mg/ml (+/-10 %)                                                                                                          |
| <b>Conjugation</b>    | Unconjugated                                                                                                                 |
| <b>Purification</b>   | Purity accurate to greater than or equal to 95%                                                                              |
| <b>Dilution Range</b> |                                                                                                                              |
| <b>Formulation</b>    | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative                                                         |
| <b>Isotype</b>        | IgG1                                                                                                                         |
| <b>Storage</b>        | 24 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21 |
| <b>Instruction</b>    | OK +35 °C, 21 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer.                       |

**TARGET INFORMATION**

|                           |                 |
|---------------------------|-----------------|
| <b>Gene ID</b>            |                 |
| <b>Gene Symbol</b>        |                 |
| <b>Uniprot ID</b>         |                 |
| <b>Immunogen</b>          | human C-peptide |
| <b>Immunogen Region</b>   |                 |
| <b>Specificity</b>        |                 |
| <b>Immunogen Sequence</b> |                 |