

Anti-Chlamydia trachomatis antibody [6709] (STJ400038)

STJ400038

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Chlamydia trachomatis is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	C.trachomatis

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6709
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG2a
Storage	36 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21
Instruction	OK +35 °C, 21 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Chlamydia trachomatis LPS
Immunogen Region	
Specificity	
Immunogen Sequence	