

Anti-CYFRA21-1 antibody [1602] (STJ400057)

STJ400057

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CYFRA21-1 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1602
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.05 % SulfoBetaine, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage	6 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21 OK
Instruction	+35 °C, 21 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	human cytokeratin 19 CYFRA21-1 fragment
Immunogen Region	
Specificity	
Immunogen Sequence	