

Anti-CG antibody [5009] (STJ400084)

STJ400084

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CG is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5009
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	18 months from manufacturing at 2°C

TARGET INFORMATION

Gene ID	1511
Gene Symbol	CTSG
Uniprot ID	CATG_HUMAN
Immunogen	human chorionic gonadotropin
Immunogen Region	
Specificity	
Immunogen Sequence	