

**Anti-HE4 antibody [4503] (STJ400096)**

STJ400096

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-HE4 is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	4503
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purity accurate to greater than or equal to 95%
<b>Dilution Range</b>	
<b>Formulation</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
<b>Isotype</b>	IgG1
<b>Storage</b>	Unspecified, storage at 2°C. Temperature analysis:-70 °C, 21 days OK-20 °C, 21 days OK +4 °C, 21 days OK +35 °C, 7 days OK +35 °C, 21 days OK +45 °C, 3 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">10406</a>
<b>Gene Symbol</b>	<a href="#">WFDC2</a>
<b>Uniprot ID</b>	<a href="#">WFDC2_HUMAN</a>
<b>Immunogen</b>	human epididymis protein 4 (HE4)
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	