

Anti-LH antibody [5302] (STJ400127)

STJ400127

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-LH is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5302
Concentration	1.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage	24 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days N/D +4 °C, 21 days N/D +35 °C, 7
Instruction	days N/D +35 °C, 21 days N/D +45 °C, 3 days N/D +45 °C, 7 days N/D. Shelf life stability testing at 2–8 °C in the product bu

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	human luteinizing hormone (lutropin) , and its beta-subunit
Immunogen Region	
Specificity	
Immunogen Sequence	