

## Anti-KLK3 antibody [8311] (STJ400167)

STJ400167

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-KLK3 is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	8311
<b>Concentration</b>	1.0-1.1mg/ml
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purity accurate to greater than or equal to 95%
<b>Dilution Range</b>	
<b>Formulation</b>	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
<b>Isotype</b>	IgG2a
<b>Storage</b>	36 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21
<b>Instruction</b>	OK +35 °C, 21 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">354</a>
<b>Gene Symbol</b>	<a href="#">KLK3</a>
<b>Uniprot ID</b>	<a href="#">KLK3_HUMAN</a>
<b>Immunogen</b>	human prostate-specific antigen
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	