

Anti-SHBG antibody [6001] (STJ400175)

STJ400175

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SHBG is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6001
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage	24 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21
Instruction	days OK +35 °C, 7 days OK +35 °C, 21 days Failed due to reduced antigen binding +45 °C, 3 days OK +45 °C, 7 days OK, some cha

TARGET INFORMATION

Gene ID	6462
Gene Symbol	SHBG
Uniprot ID	SHBG_HUMAN
Immunogen	human sex hormone-binding globulin
Immunogen Region	
Specificity	
Immunogen Sequence	