

Anti-ST2 antibody [10201] (STJ400210)

STJ400210

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-ST2 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10201
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage	Unspecified, storage at 2–8°C. Temperature analysis: -70 °C, 21 days OK -20 °C, 21 days OK +4 °C, 21 days OK +35 °C, 21 days OK
Instruction	+45 °C, 3 days OK +45 °C, 7 days Minor charge alterations. Shelf life stability testing at 2–8 °C in the product buffer.

TARGET INFORMATION

Gene ID	6822
Gene Symbol	SULT2A1
Uniprot ID	ST2A1_HUMAN
Immunogen	Human soluble ST2
Immunogen Region	
Specificity	
Immunogen Sequence	