

Anti-PINP antibody [11502] (STJ400257)

STJ400257

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-PINP is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11502
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage	Unspecified, storage at -20°C. Shipping at 2–8°C.. Temperature analysis: -70 °C, 21 days Not Determined (N/D) -20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Intact form of Human procollagen I N-terminal peptide
Immunogen Region	
Specificity	
Immunogen Sequence	