

Anti-Gastrin-17 antibody [11701] (STJ400258)

STJ400258

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Gastrin-17 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11701
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	
Storage	Temperature analysis:-70 °C, 21 days Not Determined (N/D)-20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7
Instruction	days OK. Shelf life stability testing at 2-8 °C in the product buffer.

TARGET INFORMATION

Gene ID	2520
Gene Symbol	GAST
Uniprot ID	GAST_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	