

Anti-Calcitonin antibody [4051] (STJ400271)

STJ400271

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Calcitonin is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4051
Concentration	> 1.0 mg/ml
Conjugation	Unconjugated
Purification	Greater than 95% purity
Dilution Range	
Formulation	0.1 M BBS, pH 8, 0.9 % NaCl
Isotype	IgG1
Storage Instruction	60 months from manufacturing at 2–8°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human calcitonin and procalcitonin
Immunogen Region	
Specificity	
Immunogen Sequence	