

**Anti-Renin antibody [10850] (STJ400273)**

STJ400273

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Renin is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	10850
<b>Concentration</b>	> 1.0 mg/ml
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN <sub>3</sub> as a preservative
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	60 months from manufacturing at 2–8°C.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">5972</a>
<b>Gene Symbol</b>	<a href="#">REN</a>
<b>Uniprot ID</b>	<a href="#">REN1_HUMAN</a>
<b>Immunogen</b>	Human renin
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	