

Anti-Renin antibody [10852] (STJ400275)

STJ400275

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Renin is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10852
Concentration	> 1.0 mg/ml
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	60 months from manufacturing at 2–8°C.

TARGET INFORMATION

Gene ID	5972
Gene Symbol	REN
Uniprot ID	REN1_HUMAN
Immunogen	Human renin
Immunogen Region	
Specificity	
Immunogen Sequence	