

Anti-Albumin antibody [6501] (STJ400279)

STJ400279

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Albumin is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6501
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	36 months from manufacturing at 2–8 °C

TARGET INFORMATION

Gene ID	213
Gene Symbol	ALB
Uniprot ID	ALBU_HUMAN
Immunogen	N/D
Immunogen Region	
Specificity	Antibody recognizes human albumin
Immunogen Sequence	