

Anti-CK-MB antibody [7501] (STJ400282)

STJ400282

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CK-MB is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	7501
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	24 months from manufacturing at 2–8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	N/D
Immunogen Region	
Specificity	Antibody recognizes human creatine kinase MB isoenzyme
Immunogen Sequence	