

**Anti-GDF-15 antibody [4901] (STJ400288)**

STJ400288

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-GDF-15 is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	4901
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Unspecified, storage at 2-8 °C

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">9518</a>
<b>Gene Symbol</b>	<a href="#">GDF15</a>
<b>Uniprot ID</b>	<a href="#">GDF15_HUMAN</a>
<b>Immunogen</b>	N/D
<b>Immunogen Region</b>	
<b>Specificity</b>	Antibody recognizes human growth differentiation factor 15 (GDF-15)
<b>Immunogen Sequence</b>	