

Anti-GDF-15 antibody [4903] (STJ400290)

STJ400290

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-GDF-15 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4903
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2-8 °C

TARGET INFORMATION

Gene ID	9518
Gene Symbol	GDF15
Uniprot ID	GDF15_HUMAN
Immunogen	N/D
Immunogen Region	
Specificity	Antibody recognizes human growth differentiation factor 15 (GDF-15)
Immunogen Sequence	