

Anti-NT-proBNP antibody [RC1306] (STJ400297)

STJ400297

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse, Human monoclonal antibody anti-NT-proBNP is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse, Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RC1306
Concentration	5.0 mg/ml (+/-10 %) (A280 nm, 1 mg/ml, 1 cm = 1.6)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	human IgG1
Storage Instruction	Unspecified, storage at 2–8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Amino acids 12-28 (LETSGLQEQRNHLQGKL)
Immunogen Region	
Specificity	Antibody recognizes human N-terminal pro-B-type natriuretic peptide
Immunogen Sequence	