

Anti-SAA antibody [2205] (STJ400309)

STJ400309

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SAA is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2205
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2-8 °C

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recognizes aa 64-71 but also reacts with aa 28-38 indicating that the epitope may be partly conformational.
Immunogen Region	
Specificity	Antibody recognizes human serum amyloid A (SAA)
Immunogen Sequence	