

Anti-ACTH antibody [R12802] (STJ400341)

STJ400341

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	NA monoclonal antibody anti-Adrenocorticotrophic Hormone is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R12802
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Liquid, may turn slightly opaque during storage
Isotype	IgG1
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	0
Immunogen Region	
Specificity	Antibody recognizes human adrenocorticotrophic hormone
Immunogen Sequence	