

**Anti-BNP antibody [11906] (STJ400343)**

STJ400343

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	NA monoclonal antibody anti-Brain Natriuretic Peptide is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	11906
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Liquid, may turn slightly opaque during storage
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	N/D (N-terminal or mid-region based on pairing information)
<b>Immunogen Region</b>	
<b>Specificity</b>	Antibody recognizes human B-type natriuretic peptide
<b>Immunogen Sequence</b>	