

Anti-BNP antibody [11908] (STJ400344)

STJ400344

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | NA monoclonal antibody anti-Brain Natriuretic Peptide is suitable for use in ELISA research applications. |
| Applications | ELISA |
| Host/Source | |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | 11908 |
| Concentration | 5.0 mg/ml (+/-10 %) |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | Liquid, may turn slightly opaque during storage |
| Isotype | IgG1 |
| Storage Instruction | |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | Amino acid region 10–21. |
| Immunogen Region | |
| Specificity | Antibody recognizes human B-type natriuretic peptide |
| Immunogen Sequence | |