

Anti-HBP antibody [12201] (STJ400345)

STJ400345

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	NA monoclonal antibody anti-Heparin Binding Protein is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	12201
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	