

**Anti-PIVKA-II antibody [12103] (STJ400349)**

STJ400349

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	NA monoclonal antibody anti-PIVKA-II is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	12103
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Liquid, may turn slightly opaque during storage
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Amino acid region 23-46 (SYEEAFEALSTATDVFWAKYTA)
<b>Immunogen Region</b>	
<b>Specificity</b>	Antibody recognizes human PIVKA-II (Protein Induced Vitamin K absence or Antagonist-II)
<b>Immunogen Sequence</b>	