

Anti-HTR2C antibody (400-458) {Biotin} (STJ500010)

STJ500010

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HTR2C (400-458) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:100 ICC: 1:50-1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3358
Gene Symbol	HTR2C
Uniprot ID	5HT2C_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-458 on human 5-hydroxytryptamine receptor 2C.
Immunogen Region	400-458
Specificity	
Immunogen Sequence	