

Anti-ANGPTL4 antibody (50-100) {Biotin} (STJ500082)
STJ500082**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ANGPTL4 (50-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57875
Gene Symbol	Angptl4
Uniprot ID	ANGL4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 50-100 on mouse Angiopoietin-related protein 4.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	