

Anti-ANGPTL5 antibody (340-388) (STJ500084)

STJ500084

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ANGPTL5 (340-388) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.66 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	253935
Gene Symbol	ANGPTL5
Uniprot ID	ANGL5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 340-388 on human ANGPTL5 protein.
Immunogen Region	340-388
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat small intestine with Anti-ANGPTL5 primary antibody (STJ500084) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.