

**Anti-AQP11 antibody (140-190) {FITC} (STJ500122)**  
STJ500122**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-AQP11 (140-190) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM/ELISA/ICC/IF/IHC/IP/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Pig/Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.57-0.67 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution</b>	WB: 1:500
<b>Range</b>	ELISA: 1:10, 000
	IP: 1:200
	IHC: 1:500
	ICC: 1:50
<b>Formulation</b>	IF: 1:50
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">286758</a>
<b>Gene Symbol</b>	<a href="#">Aqp11</a>
<b>Uniprot ID</b>	<a href="#">AQP11_RAT</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 140-190 on Aquaporin 11 protein.
<b>Immunogen Region</b>	140-190
<b>Specificity</b>	
<b>Immunogen Sequence</b>	