

Anti-Baboon antibody (552-572) (STJ500202)

STJ500202

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Baboon (552-572) is suitable for use in Confocal Microscopy, ELISA, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM/ELISA/IF/IHC/IP/WB
Host/Source	Rabbit
Reactivity	D.melanogaster

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.66 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:10, 000
	ELISA: 1:10, 000
	IP: 1:200
	IHC: 1:200
	ICC: 1:200
	IF: 1:200
Formulation	CFM: 1:200
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	C-epitope Synthetic peptide corresponding to amino acids 552-572 of Drosophila melanogaster Baboon protein. Peptide was covalently modified.
Immunogen Region	552-572
Specificity	
Immunogen Sequence	