

**Anti-BTBD7 antibody (1050-1150) {FITC} (STJ500302)**

STJ500302

**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-BTBD7 (1050-1150) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications. |
| <b>Applications</b>      | CM/ELISA/ICC/IF/IHC/IP/WB   |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Human/Mouse   |

**PRODUCT PROPERTIES**

|                      |  |
|----------------------|--|
| <b>Clonality</b>     | Polyclonal   |
| <b>Clone ID</b>      |  |
| <b>Concentration</b> | 0.61-0.68 µg/µl  |
| <b>Conjugation</b>   | FITC   |
| <b>Purification</b>  | Affinity Purified  |
| <b>Dilution</b>      | WB: 1:500  |
| <b>Range</b>         | DB: 1:10, 000  |
|                      | ELISA: 1:10, 000   |
|                      | IP: 1:200  |
|                      | IHC: 1:200   |
|                      | ICC: 1:200   |
| <b>Formulation</b>   | CFM: 1:200   |
|                      | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
|                      |  |
| <b>Isotype</b>       | IgG  |
| <b>Storage</b>       | Store at -20°C for long term storage. Avoid freeze-thaw cycles.  |
| <b>Instruction</b>   |  |

**TARGET INFORMATION**

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | 55727  |
| <b>Gene Symbol</b>        | BTBD7  |
| <b>Uniprot ID</b>         | BTBD7_HUMAN  |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 1050-1150 on human BTBD7 protein. |
| <b>Immunogen Region</b>   | 1050-1150  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |