

Anti-C2orf88 antibody (25-100) {Biotin} (STJ500322)

STJ500322

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C2orf88 (25-100) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.7 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:15, 000-1:20, 000 ELISA: 1:15, 000-1:20, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84281
Gene Symbol	C2orf88
Uniprot ID	SMAKA_HUMAN
Immunogen	Mid Region epitope synthetic peptide taken within amino acid region 25-100 on human C2orf88 protein.
Immunogen Region	25-100
Specificity	
Immunogen Sequence	