

Anti-CALB2 antibody (220-271) {FITC} (STJ500351)

STJ500351

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CALB2 (220-271) is suitable for use in Confocal Microscopy, Dot Blot, ELISA, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM/DB/ELISA/IF/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:500
Range	ELISA: 1:10, 000
	IP: 1:200
	IHC: 1:50-1:200
	ICC: 1:50-1:200
Formulation	IF: 1:50-1:200
	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	794
Gene Symbol	CALB2
Uniprot ID	CALB2_HUMAN
Immunogen	Synthetic peptide within amino acid region 220-271. Conjugated to Keyhole Limpet Haemocyanin (KLH). Peptide was covalently modified.
Immunogen	220-271
Region	
Specificity	
Immunogen	
Sequence	