

Anti-CCR4 antibody (1-50) {Biotin} (STJ500406)

STJ500406

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CCR4 (1-50) is suitable for use in Confocal Microscopy, Dot Blot, ELISA, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CF/DB/ELISA/IF/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Cat/Human/Monkey/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution	WB: 1:500
Range	IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1233
Gene Symbol	CCR4
Uniprot ID	CCR4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human C-C chemokine receptor type 4 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	