

Anti-CXCR2 antibody (18-39) {Biotin} (STJ500650)

STJ500650

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CXCR2 (18-39) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Flow Cytometry and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/IP/FC/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.60 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:250 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3579
Gene Symbol	CXCR2
Uniprot ID	CXCR2_HUMAN
Immunogen	Synthetic peptide: EDLSNYSYSTLPFLDAAAPC Corresponding to amino acids 18-39 of Human CXCR2.
Immunogen Region	18-39
Specificity	
Immunogen Sequence	