

Anti-Dopamine Receptor D1 antibody (290-340) {FITC} (STJ500707)

STJ500707

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dopamine Receptor D1 (290-340) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	CM/ELISA/ICC/IF/IHC/WB
Host/Source	Rabbit
Reactivity	Human/Monkey/Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-1.20 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	Drd1
Uniprot ID	DRD1_RAT
Immunogen	Synthetic peptide taken within amino acid region 290-340 on human D-1 dopamine receptor protein.
Immunogen Region	290-340
Specificity	
Immunogen Sequence	